

Gebrauchsanweisung / Instructions for use / Istruzioni per l'uso / Mode d'emploi / Gebruiksaanwijzing / Instrucciones de uso / Käyttöohjeet / Bruksanvisning / Bruksanvisning / Инструкции за употреба / Οδηγίες χρήσης / Uputstvo za upotrebu / Návod k použití / Instrukcja obsługi / Інструкція із застосування / Instrucțiuni de utilizare / Návod na použitie / Navodila za uporabo / Használati utasítás / Upute za korištenje / Naudojimo instrukcijos / Lietošanas instrukcija / Kasutusjuhend / Инструкция по применению / Kullanım kılavuzu / Instruções de uso /取扱説明書 / 사용 설명서

DE	XD 2035	Patientenschlauch 1,50m	/ XD 2040	Patientenschlauch 2,50m	/ XD 2045	Patientenschlauch 3,20m
EN	XD 2035	Patient tubing 1.50 m	/ XD 2040	Patient tubing 2.50 m	/ XD 2045	Patient tubing 3.20 m
IT	XD 2035	Tubo per paziente 1,50 m	/ XD 2040	Tubo per paziente 2,50 m	/ XD 2045	Tubo per paziente 3,20 m
FR	XD 2035	Tubulure patient 1,50 m	/ XD 2040	Tubulure patient 2,50 m	/ XD 2045	Tubulure patient 3,20 m
NL	XD 2035	patiëntslang 1,50 m	/ XD 2040	patiëntslang 2,50 m	/ XD 2045	patiëntslang 3,20 m
ES	XD 2035	Tubo del paciente 1,50 m	/ XD 2040	Tubo del paciente 2,50 m	/ XD 2045	Tubo del paciente 3,20 m
FI	XD 2035	potilasletku 1,50 m	/ XD 2040	potilasletku 2,50 m	/ XD 2045	potilasletku 3,20 m
NO	XD 2035	pasientslange, 1,50 m	/ XD 2040	pasientslange, 2,50 m	/ XD 2045	pasientslange, 3,20 m
SV	XD 2035	Patientslang 1,50 m	/ XD 2040	Patientslang 2,50 m	/ XD 2045	Patientslang 3,20 m
BG	XD 2035	Пациентен маркуч 1,5 m	/ XD 2040	Пациентен маркуч 2,5 m	/ XD 2045	Пациентен маркуч 3,2 m
EL	XD 2035	Καθετήρας 1,50 m	/ XD 2040	Καθετήρας 2,50 m	/ XD 2045	Καθετήρας 3,20 m
SR	XD 2035	Crevo za pacijenta 1,50 m	/ XD 2040	Crevo za pacijenta 2,50 m	/ XD 2045	Crevo za pacijenta 3,20 m
CS	XD 2035	pacientská hadička 1,50 m	/ XD 2040	pacientská hadička, 2,50 m	/ XD 2045	pacientská hadička, 3,20 m
PL	XD 2035	Wężyki pacjenta 1,50 m	/ XD 2040	Wężyki pacjenta 2,50 m	/ XD 2045	Wężyki pacjenta 3,20 m
UA	XD 2035	Система трубок пацієнта 1,5 м	/ XD 2040	Система трубок пацієнта 2,5 м	/ XD 2045	Система трубок пацієнта 3,2 м
RO	XD 2035	Tubulatură pentru pacient 1,50 m	/ XD 2040	Tubulatură pentru pacient 2,50 m	/ XD 2045	Tubulatură pentru pacient 3,20 m
SK	XD 2035	Hadička pacienta 1,50 m	/ XD 2040	Hadička pacienta 2,50 m	/ XD 2045	Hadička pacienta 3,20 m
SL	XD 2035	Cevke za bolnika 1.50 m	/ XD 2040	Cevke za bolnika 2,50 m	/ XD 2045	Cevke za bolnika 3,20 m
HU	XD 2035	páciensszerezélék 1,50 m	/ XD 2040	páciensszerezélék 2,50 m	/ XD 2045	páciensszerezélék 3,20 m
HR	XD 2035	Cijev za bolesnika 1,50 m	/ XD 2040	Cijev za bolesnika 2,50 m	/ XD 2045	Cijev za bolesnika 3,20 m
LT	XD 2035	paciento vamzdelis 1,50 m	/ XD 2040	paciento vamzdelis 2,50 m	/ XD 2045	paciento vamzdelis 3,20 m
LV	XD 2035	Pacienta caurulītes 1,50 m	/ XD 2040	Pacienta caurulītes 2,50 m	/ XD 2045	Pacienta caurulītes 3,20 m
ET	XD 2035	patsiendivoolik 1,50 m	/ XD 2040	patsiendivoolik 2,50 m	/ XD 2045	patsiendivoolik 3,20 m
RU	XD 2035	Трубка пациента, длина 1,50 м (XD 2035)	/ XD 2040	Трубка пациента, длина 2,50 м (XD 2040)	/ XD 2045	Трубка пациента, длина 3,20 м (XD 2045)
TR	XD 2035	Hasta hortumu 1,50 m	/ XD 2040	Hasta hortumu 2,50 m	/ XD 2045	Hasta hortumu 3,20 m
BR	XD 2035	Tubo flexível de paciente 1,50 m	/ XD 2040	Tubo flexível de paciente 2,50 m	/ XD 2045	Tubo flexível de paciente 3,20 m
JA	XD 2035	患者チュービング 1.50 m	/ XD 2040	患者チュービング 2.50 m	/ XD 2045	患者チュービング 3.20 m
KO	XD 2035	환자 튜브 1.50m	/ XD 2040	환자 튜브 2.50m	/ XD 2045	환자 튜브 3.20m

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Abkürzungen

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EN Instructions for use

Abbreviations

Abbreviation	Significance
CT	Computed tomography
MRI	Magnetic Resonance Imaging
CM	Contrast medium
NaCl	Saline solution (0.9% saline solution)
PET	Positron Emission Tomography

About This Document

- These instructions for use are part of the product and they describe its safe intended use.
- Read the instructions for use before using the product.
- Keep the instructions for use with the product so that they are accessible at all times.

Labeling and symbols

Marking	Significance
 WARNING	Warning! Hazards to persons. Non-compliance can cause death or serious injuries to persons, or damage to the product.
→	Handling instructions

Symbols on the product and packaging

Icon	Description	Icon	Description
	Consult instructions for use		Sterilized using ethylene oxide
	CE mark with identification number of the notified body		Do not use if package is damaged and consult instructions for use
	Catalogue number		Do not use if protective caps are not properly attached

Icon	Description	Icon	Description
	Batch code		Non-pyrogenic
	Quantity		Keep dry
	Use-by date		Temperature limit
	Manufacturer		Humidity limitation, non-condensing
	Date of manufacture		Keep away from sunlight
	Do not re-use		Caution
	Single sterile barrier system		Medical device
	Patient tubing		

Intended Purpose

Tubes used for intravenous application of liquids and drugs.

Indications for use

MRI and PET/MRI, CT and PET/CT examinations with contrast media.

Contraindications

Administration of contrast media during high-pressure angiography. Other applications that do not comply with the intended use of the patient hose.

Side effects

No side effects are known in association with the pump tubing. The residual risks defined in the risk management are covered by the safety notes in the instructions for use.

User profile

The patient hose may only be used by qualified medical personnel (i.e. medical technical assistants, radiologic technologists, medical assistants and radiologists) and after instruction according to the instructions for use. Instruction is provided once by ulrich medical or a person authorised by ulrich medical as part of the initial product training or product updates requiring an instruction. Depending on the hospital's internal and national regulations, authorised users may instruct other users of the same healthcare facility.

Patient group

No restrictions. Patients who receive a contrast-enhanced computer tomography or magnetic resonance imaging.

Function of medical device/Operating principle

The patient hose for ulrich medical CT/MRI injectors establishes the connection from the pump tubing to the venous access on the patient side. The media (NaCl and CM) are delivered to the patient hose in a controlled manner by a roll pump using the pump tubing.

Product description

The patient tubing consists of a tube with two protective caps (2) and has a male Luer lock connector (3) with an integrated check valve on the patient side (B). The pump tubing is connected to the female Luer lock connector (1). A second check valve (4) is integrated into the tubing. The check valves are self-closing. They open only when a medium flows in the defined direction. In case of counter-pressure or no media flow, the valves remain closed (see fig.).

- A Device side with pump tubing
- B Patient side
- 1 Female Luer lock
- 2 Protective caps for tubing ends
- 3 Male Luer lock, with the first check valve
- 4 Second check valve
- 5 Media flow direction

 **WARNING! Risk of injury due to infection, contamination, or material fatigue!**

Product intended for single use. The product is EO sterilized and delivered in sterile packaging.

- Handle all sterile components carefully to ensure sterility.
- Do not resterilize or re-use the product.
- Before use, visually check the packaging for damage and the sterile barrier for physical integrity. Do not use the product if the packaging or the sterile barrier is damaged.
- Do not use the product if the protective caps are missing or are not firmly seated.
- Do not use the product after the expiry date.
- Do not use product that was not stored under the recommended storage conditions.
- Do not use the product if any parts are damaged.
- Replace the patient tubing for each new patient.
- After disconnecting the patient tubing, the system must immediately be equipped with new patient tubing, and the protective cap must not be removed until it is connected to the patient.
- Dispose of patient tubing upon disconnecting it from the patient in order to prevent contamination of the pump tubing.

Safe handling

- ➔ Before using the product, check that it is operable and in proper condition. Do not use the product if any parts are damaged.
- ➔ In order to avoid damage due to improper use, use the product in accordance with these instructions for use and the injector's instructions for use.
- ➔ Product and accessories must only be used by persons who have the required knowledge, training or experience.
- ➔ Ensure that all connectors are correctly connected to each other. DO NOT OVERTIGHTEN, because this can break the threads of the Luer lock connector and lead to ruptures.
- ➔ Being closed too tightly and being in prolonged contact with contrast media can make the Luer lock connection hard to unscrew. Carefully unscrew the Luer lock connector.

WARNING! Risk of injury! Flow reduction through pressure regulation, air injection and overdosage through the use of tubing extensions!

The use of tubing extensions and tubing combinations or systems that have not been approved by ulrich medical jeopardizes the safety of the patient and risks the contamination of the tubing system.

- ➔ Do not use tubing extensions (e.g. "Heidelberg extension") and tubing combinations or systems that have not been approved by ulrich medical.
- ➔ Only use original products approved by ulrich medical.

Reporting incidents

All complaints in connection with the safety, efficacy or performance of the product must be reported by the user to ulrich medical (complaint@ulrichmedical.com) or the local distributor within 72 hours.

If one or more components of the system should have a malfunction (i.e. one or more of the performance specifications are not met or the intended performance is not achieved for other reasons) or if it is suspected that this may be the case, (complaint@ulrichmedical.com) or the local distributor must be notified immediately (within 24 hours).

If one or more components of the system has ever failed and possibly led or contributed to the death or severe injury of a patient, user or third party, then ulrich medical or the local distributor must be notified by telephone immediately. Please use the Medical Device Vigilance Report Form from ulrich medical to report product complaints, malfunctions or product failure. It is available at: www.ulrichmedical.de/meldebogen.

Use

- ➔ Remove the product from the packaging and check for damaged parts before the product is used.
- ➔ Connect the product to the pump tubing and use it according to the instructions in the injector user manual.
- ➔ Only use the product with ulrich medical original products*:

Art. No.	Designation
XD 8000	ulrichINJECT CT motion
XD 8151	Pump tubing
XD 8003	Pump tubing

* The country-specific availability of items must be taken into account.

Technical data

Parameter	Value
Pressure tightness	Max. 22.4 bar
Material	PVC (DEHP free, latex free)
Max. service life	12 hours for one patient

Storage

- ➔ Store the product in the sterile packaging in a dry, dark and temperature-controlled room.
- ➔ Protect the sterile packaged product from environmental factors such as dirt, dust, humidity, heat, and sunlight.

Ambient conditions

	Use	Storage
Temperature	+15°C... +40°C	+15°C... +30°C
Relative humidity	10% ... 90% non-condensing	30% ... 60% non-condensing

Disposal

- ➔ After use, dispose of the product and packaging in accordance with the local applicable regulations for disposable products contaminated with blood. Of course, the system and system components can also be disposed of via the clinic's own disposal system. The disposal of medical waste in the European Union is subject to the European List of Waste Materials. If necessary, various levels of government enact their own waste regulations, which must be observed. Outside the European Union, the relevant country-specific regulations for the disposal of medical waste must be observed.

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Ulrich GmbH & Co. KG | Buchbrunnenweg 12 | 89081 Ulm | Germany